

Understanding variability in microdialysis measurements: introducing a combined calibration approach for piperacillin and tazobactam in LPS-induced septic piglets

F. Mueller (1,2), E. Hermans (3,4,5), D. Bindellini (1,2), R. Michelet (1), L.B.S. Aulin (1), P. De Paepe (3,6), P. De Cock (3,7,8), M. Devreese (4), C. Klotz (1,2)

(1) Department of Clinical Pharmacy and Biochemistry, Institute of Pharmacy, Freie Universität Berlin, Germany, (2) and Graduate Research Training program PharMetX, Germany, (3) Department of Basic and Applied Medical Sciences, Ghent University, Belgium, (4) Department of Pathobiology, Pharmacology and Zoological Medicine, Ghent University, Belgium, (5) Department of Pediatrics, Ghent University Hospital, Belgium, (6) Department of Emergency Medicine, Ghent University Hospital, Belgium, (7) Department of Pharmacy, Ghent University Hospital, Belgium, (8) Department of Pediatric Intensive Care, Ghent University Hospital, Belgium

Background and Objectives

Inadequate target tissue fluid exposure may lead to treatment failure [1] in septic children → impact of **sepsis on tissue fluid exposure** in children for **piperacillin (PIP)** and **tazobactam (TAZ)** not investigated

Lipopolysaccharide (LPS) induced **septic piglet model** can elucidate impact of septic like state on drug disposition and pharmacokinetics

Microdialysis (μD): Gold standard method to obtain unbound drug **interstitial space fluid concentrations** (C_{ISF}) [2], **reliable catheter calibration** is required for meaningful results

Catheter-specific calibration is commonly conducted by (i) **retrodialysis** (before/after dosing drug of interest) or (ii) an **internal standard** parallel to microdialysis sampling (**IS**), with penicillin (PEN) used for PIP and TAZ [3]

NLME methods are the **most adequate approach** to handle **μD-based data analysis** [4], IS calibration has never been included yet

We investigated the **feasibility** and potential **benefits** of both calibration methods individually and their combination leveraging μD-data from a **porcine animal study**

Methods

Calibration methods

- Retrodialysis (rD), Fig. 1
- Internal standard (IS)

Catheter insertion

- 1 catheter per piglet
- Into paraspinal musculature ISF

Study design

- Dosing (5 per period, tau 6 h)
- PIP 75 mg/kg; TAZ 9.375 mg/kg
- LPS infusion on day 3
- Sampling after 1st and 5th dose (steady-state) per period




Figure 1: Principle of microdialysis and retrodialysis technique; exemplified by interstitial space fluid in adipose tissue.




Figure 2: Study design

Results

Plasma model

- Separate** models for PIP and TAZ (Figure 4)
- Retrodialysis as reference model
- Covariates** to account for **septic-like state** included (Figure 4)
- Plasma-derived **PK parameters** (Table 1) **remained constant**: Change <10% compared to retrodialysis model
- RSE below <65%

Parameter [unit]	PIP estimates	TAZ estimates
CL [L/h]	4.01	5.90
V1 [L]	1.35	1.88
Q [L/h]	12.5	18.3
V2 [L]	1.62	1.80
IIV and IOV, %CV		
IIV on CL	10.0	8.38
IIV on V1	24.2	18.0
IOV on CL	15.7	9.40
Residual variability		
add C_{Plasma} [mg/L]	0.56	0.083
prop C_{Plasma} [CV.%]	16.5	17.0

Retrodialysis

Parameter [unit]	PIP estimates	TAZ estimates
rD [%]	23.2	29.4
TF [%]	115	139
Time_rD	0.044	0.039
rD _{PEN} _rD	-	-
IIV and IOV, %CV		
IIV on rD	14.5	12.0
IOV on rD	16.8	13.2
Residual variability		
prop μD [CV.%]	22.7	17.8
add μD [mg/L]	0.27	0.29
prop CrD [CV.%]	10.8	12.0

IS calibration

PIP estimates	TAZ estimates
29.6*	29.5*
88.3	136
0.0013	0.0011
-	-
10.4	8.5
20.3	19.4
25.8	22.4
0.24	0.30
5.45	10.7

Combined approach

PIP estimates	TAZ estimates
22.6	28.7
117	140
0.0241	0.026
1.26	1.19
8.7	6.6
16.0	12.9
21.9	17.0
0.271	0.28
10.1	11.3

Calibration

- Increased rD in control group over study days
- Time effect on rD for control group: $rD_{TIME} = (((Time)/40) * \theta(Time_rD))$
- Decreased rD in LPS septic like state (Figure 5)




Figure 5: rD values over study days for piperacillin, tazobactam and penicillin per group with their medians

Discussion and Conclusions

Retrodialysis calibration not considering physiological changes while sampling
→ Assumption: Catheter extracts same fraction throughout entire sampling period

IS calibration: Drug-specific differences in rD compensated by TF
No improvement for PIP/TAZ in VPC or IIV
→ Change in tissue exposure expected

Benefits combined calibration approach

- Individual drug-specific rD** for each sampling interval
- Variability** between samples **explained by IS**
- Improved **predictive performance** in control group
- Reduction of rD-related IIV**
- Reduction rD-related time dependency**
- Feasible and beneficial**

Changing **blood flow** though pigs' **activity** could explain
1) Increase of general rD over days
2) Variability in rD_{PEN} over time
→ **Blood flow** and **activity** should be **monitored**

IS calibration in combination with retrodialysis explained *in vivo* observed rD-related variability, thus improves understanding and interpretation of μD data

References

- C. Fleischmann-Struzek et. al, Lancet Respir Med, 2018
- N. Plock et. al, Eur. J. Pharm. Sci., 2005
- M.B. Knudsen et. al, J. Pharm. Sci. 2021
- D. Busse et. al, Pharm.Res., 2021

Abbreviations

Add: additive; ALB: Albumin; CL: Clearance; CRT: Creatinine IIV: Interindividual variability; IOV: Inter-occasional variability; ISF: Interstitial space fluid; LEU: Leucocytes; LPS: Lipopolysaccharides; PIP: Piperacillin; PEN: Penicillin; prop: proportional RR: Relative recovery; rD: Relative delivery; TAZ: Tazobactam; PEN: Penicillin; TEMP: Temperature; Time_rD: Time effect on rD; rD_{PEN}_rD: fractional change of rD based on IS calibration; V: Volume of distribution; VPC: Visual predictive check; WGT: Body weight; μD: Microdialysis



For more information:
Felix.leon.mueller@fu-berlin.de
www.clinical-pharmacy.eu

32nd Population Approach Group Europe-meeting - PAGE, Rome, Italy 2024

Poster PDF



AK Klotz

